


PRESSURE MYOGRAPH SYSTEM – 110P-XL



- MULTI-PURPOSE FUNCTIONALITY
- COMPACT DESIGN
- RELIABLE & ROBUST MATERIALS
- ACCURATE & RELIABLE DATA
- SMALL FOOT-PRINT & PORTABLE
- DIRECT DATA STREAMING
- STATE-OF-THE-ART TECHNOLOGY

The Pressure Myograph System – 110P-XL is used to study the structure and function of large arteries, veins, and other vessels. Large vessel reactivity and compliance can be measured during diameter changes while the tissue is under physiological pressures, precisely controlled via the system interface.

The tissue mounts can accommodate vessels with internal diameters of 1.5 – 6.0 mm. A built-in heating system maintains the chamber temperature, eliminating the need for continuous and often costly superfusion. The chamber cover includes ports for superfusion, rapid draining and filling, oxygenation and cumulative addition of drugs. To facilitate cleaning, the chamber is made of acid resistant stainless steel.

Once the vessel is mounted onto the two cannulas, the large optical window allows for tissue observation and continuous diameter measurements. The dimension analysis software MyoVIEW will control the temperature and internal pressure, while also recording diameter changes.

The 110P-XL system generates constant pressures, so pressures can be adjusted, but flow is not possible. Therefore the primary area of research is physiological responses under myogenic tone. Any experiment studying a pressurized vessel can be performed.

The state of the art Pressure Myograph software MyoVIEW will collect data such as changes in outer wall diameter, intravascular pressure, and a host of other calculated parameters, setting pressure or perfusion myography apart from standard organ bath techniques. However, the use of low-power objectives will have an impact on the ability to image through the vessel to the internal walls.

The Acquisition & Analysis Package
The DMT Inverted Microscope, inverted Zeiss, Nikon or similar inverted microscopes (contact DMT for further specifications) with USB camera, computer and Data Acquisition Software - MyoVIEW.



SYSTEM SPECIFICATIONS

CHAMBER

Chamber Volume (min)	8.0 ml
Chamber Volume (max)	20.0 ml
Chamber	1
Chamber Material	Stainless steel – acid resistant
Vessel Size (diameter)	1.5 mm- 6.0 mm
Vessel Size (length)	3.0 mm- 33.0 mm
Vessel alignment	X, Y, Z
Micrometer resolution	0.01 mm
Mounting type	Cannulas (1.5mm, 2.0mm, 3.0mm and 5.0mm)

TEMPERATURE

Temperature Range	15.0 to 50.0 °C
Temperature Resolution	0.1 °C
Temperature Stability	±0.2 °C
Heating	YES

TRANSDUCER PRESSURE

Pressure Output	mmHg
Pressure Range	0 - 250 mmHg
Pressure stability	±0.5 mmHg
Pressure Resolution	0.1 mmHg

RESERVOIR

Reservoir Heated	YES
Reservoir capacity	250ml
Pressure circuit	Closed
Pressure Circuit	1 bar (max)

OUTPUT

Data communication	USB 2.0
Analogue channels	4
Analogue output range	±2.5 V

ADD-ONS AVAILABLE

- MyoVIEW Acquisition
- Computer Package (Desktop)
- Inverted Microscope
- Stage Micrometer
- Camera – High-Speed
- Vacuum Package
- Heating Bath Circulator
- Dissection Stereo Microscope
- Light Source - Arm Swan-neck
- Dissection Tool Kit
- Nylon Suture
- Gas supply manifold- 4-channel
- Peristaltic pump
- Cannulas – Other sizes optional
- Service & Maintenance
- Agreement
- Extended warranty +3 years *

* 2-Years of Warranty included in purchase